




CE0197

## Disposable A.V. Fistula Needle Sets



Gauge	Needle Length	Wing	Tube Spec	Sterilization
14G	25±2mm	fixed	ID3.2-OD5.0	EO/Gamma
14G	25±2mm	rotatable	ID3.2-OD5.0	EO/Gamma
14G	32±2mm	fixed	ID3.2-OD5.0	EO/Gamma
14G	32±2mm	rotatable	ID3.2-OD5.0	EO/Gamma
15G	25±2mm	fixed	ID3.2-OD5.0	EO/Gamma
15G	25±2mm	rotatable	ID3.2-OD5.0	EO/Gamma
15G	32±2mm	fixed	ID3.2-OD5.0	EO/Gamma
15G	32±2mm	rotatable	ID3.2-OD5.0	EO/Gamma
16G	25±2mm	fixed	ID3.2-OD5.0	EO/Gamma
16G	25±2mm	rotatable	ID3.2-OD5.0	EO/Gamma
16G	32±2mm	fixed	ID3.2-OD5.0	EO/Gamma
16G	32±2mm	rotatable	ID3.2-OD5.0	EO/Gamma
17G	25±2mm	fixed	ID3.2-OD5.0	EO/Gamma
17G	25±2mm	rotatable	ID3.2-OD5.0	EO/Gamma
17G	32±2mm	fixed	ID3.2-OD5.0	EO/Gamma
17G	32±2mm	rotatable	ID3.2-OD5.0	EO/Gamma

- Ultra thin and siliconized wall for smooth penetration
- Special process on blade reduces pain and skin damage
- Medical grade PVC tubing minimizes side effect to health

### BAIN MEDICAL EQUIPMENT(GUANGZHOU)CO.,LTD.

**Add:** No.10 Juncheng Road ,Eastern Area,Economic and Technological Development District,  
Guangzhou,510760 China

**Tel:** +86-(0)20-66856868

**Fax:** +86-(0)20-32067500

**E-mail:** sales@baingz.com

www.bainmedical.com